

pH in soil.

1. Principle

This method is an electrometric procedure for measuring pH in soils based on EPA METHOD 9045D Soil and waste pH.

2. Sample preparation

Take the soil sample...

3. Standard Preparation and Calibration

Calibrate the pH meter according to the instruments manual.

4. Method

To 20 gram of soil in a 50 ml centrifuge tube, add 20 ml of dem. water close and shake/mix it for 5 minutes.

Centrifuge it 12 minutes at 4200 RPM.

Additional dilutions are allow if working with hygroscopic soils and salts or other problematic matrices.

Measure the pH in the supernatant.

5. Equipment

pH meter with a glass electrode

Centrifuge

Mixer

Scale

Centrifuge tube (50 ml)

Equipment for taking the soil sample

6. Chemicals

pH buffers the pH meter require.

7. Safety

Always work as the principles in **Good Laboratory Practice** prescribe.